









Parameter	Preservative	Sample Holding Time	Suggested Sample Size	Type of Container	Sampling Instructions
Benzo(a)pyrene	Cool 6°C, (0.100 g) Sodium Thiosulfate 3 drops of (1:1) HCl pH<2	Sample 7 days  Extract ≤ -10°C 30 days	1 Liter	Glass with Teflon Septa 	<b>Benzo(a)pyrene Sampling Instructions</b>  7 Day Holding Time Overnight Delivery Shipped Cold <6°C  Use the cold-water faucet. Remove any faucet attachments and the aerator prior to sampling. Also, avoid rubber hose attachments.  Allow the water to run for 3 to 5 minutes prior to taking the sample to flush the water lines. Decrease the water flow to the diameter of a pencil to reduce splashing.  Do not rinse the sample bottle prior to sampling; it contains a chemical preservative. Fill the sample 1/3 full. <b>Caution: Hydrochloric Acid (HCl) is corrosive.</b> Carefully add three drops of the 1:1 Hydrochloric acid. Fill the sample bottle completely full (zero headspace) without overflowing and flushing out the preservatives. Be sure the cap is tightened to prevent leakage during shipment to the laboratory.  Sample must be collected in duplicate.  <b>RECORD BOTTLE NUMBER ON FORM</b>
Carbamates	Cool 6°C, (0.013 g) Sodium Thiosulfate, (0.370 g) Potassium Dihydrogen Citrate Buffer pH<4	28 days	40 mL	Glass with Teflon Septa 	<b>Carbamates Sampling Instructions</b>  28 Day Holding Time Overnight Delivery Shipped Cold <6°C  Use the cold-water faucet. Remove any faucet attachments and the aerator prior to sampling. Also, avoid rubber hose attachments.  Allow the water to run for 3 to 5 minutes prior to taking the sample to flush the water lines. Decrease the water flow to the diameter of a pencil to reduce splashing.  Do not rinse the sample bottle prior to sampling; it contains a chemical preservative. Fill the sample vial completely full (zero headspace) without overflowing and flushing out the preservative. Be sure the cap is tightened to prevent leakage during shipment to the laboratory.  <b>RECORD BOTTLE NUMBER ON FORM</b>
Ethylene Dibromide (EDB) Dibromochloropropane (DBCP)	Cool 6°C, (0.003 mg) Sodium Thiosulfate	Sample 14 days  Extract ≤ 4°C 24 hours	40 mL	Glass with Teflon Septa 	<b>EDB &amp; DBCP Sampling Instructions</b>  14 Day Holding Time Overnight Delivery Shipped Cold <6°C  Use the cold-water faucet. Remove any faucet attachments and the aerator prior to sampling. Also, avoid rubber hose attachments.  Allow the water to run for 3 to 5 minutes prior to taking the sample to flush the water lines. Decrease the water flow to the diameter of a pencil to reduce splashing.  Do not rinse the sample bottle prior to sampling; it contains a chemical preservative. Fill the sample vial completely full (zero headspace) without overflowing and flushing out the preservative. Be sure the cap is tightened to prevent leakage during shipment to the laboratory.  Sample must be collected in duplicate.  <b>RECORD BOTTLE NUMBER ON FORM</b>
Glyphosate	Cool 6°C, (0.004 g) Sodium Thiosulfate	14 days	40 mL	Glass with Teflon Septa 	<b>Glyphosate Sampling Instructions</b>  14 Day Holding Time Overnight Delivery Shipped Cold <6°C  Use the cold-water faucet. Remove any faucet attachments and the aerator prior to sampling. Also, avoid rubber hose attachments.  Allow the water to run for 3 to 5 minutes prior to taking the sample to flush the water lines. Decrease the water flow to the diameter of a pencil to reduce splashing.  Do not rinse the sample bottle prior to sampling; it contains a chemical preservative. Fill the sample vial completely full (zero headspace) without overflowing and flushing out the preservative. Be sure the cap is tightened to prevent leakage during shipment to the laboratory.  <b>RECORD BOTTLE NUMBER ON FORM</b>
(HAA5) Haloacetic Acid Diquat	Cool 6°C, (0.006 g) Ammonium Chloride	Sample 14 days  Extract ≤ -10°C 21 days	60 mL	Amber Glass with Teflon Septa 	<b>Haloacetic Acid Sampling Instructions</b>  14 Day Holding Time Overnight Delivery Shipped Cold <6°C  Use the cold-water faucet. Remove any faucet attachments and the aerator prior to sampling. Also, avoid rubber hose attachments.  Allow the water to run for 3 to 5 minutes prior to taking the sample to flush the water lines. Decrease the water flow to the diameter of a pencil to reduce splashing.  Do not rinse the sample bottle prior to sampling; it contains a chemical preservative. Fill the sample vial completely full (zero headspace) without overflowing and flushing out the preservative. Be sure the cap is tightened to prevent leakage during shipment to the laboratory.  Sample must be collected in duplicate.  <b>RECORD BOTTLE NUMBER ON FORM</b>
505 Package (A) Chlordane Toxaphene  505 Package (B) Aroclor(s)  505 Package (C) Organohalide(s)	Cool 6°C, (0.003 mg) Sodium Thiosulfate	Sample 14 days 7 days Heptachlor  Extract ≤ 4°C 24 hours	40 mL	Glass with Teflon Septa 	<b>505 Packages A, B, &amp; C Sampling Instructions</b>  7 Day Holding Time Overnight Delivery Shipped Cold <6°C  Use the cold-water faucet. Remove any faucet attachments and the aerator prior to sampling. Also, avoid rubber hose attachments.  Allow the water to run for 3 to 5 minutes prior to taking the sample to flush the water lines. Decrease the water flow to the diameter of a pencil to reduce splashing.  Do not rinse the sample bottle prior to sampling; it contains a chemical preservative. Fill the sample vial completely full (zero headspace) without overflowing and flushing out the preservative. Be sure the cap is tightened to prevent leakage during shipment to the laboratory.  Sample must be collected in duplicate.  <b>RECORD BOTTLE NUMBER ON FORM</b>

Parameter	Preservative	Sample Holding Time	Suggested Sample Size	Type of Container	Sampling Instructions
(VOC) Volatile Organic Compounds (TTHM) Total Trihalomethanes	Cool 6°C, (0.003 g) Sodium Thiosulfate 3 drops of (1:1) HCl pH<2  (Acid added at the time of sample collection)	14 days	40 mL	Glass with Teflon Septa 	<p><b>VOC &amp; TTHM Sampling Instructions</b></p> <p>14 Day Holding Time Overnight Delivery Shipped Cold &lt;6°C</p> <p>Use the cold-water faucet. <b>Remove any faucet attachments and the aerator prior to sampling. Also, avoid rubber hose attachments.</b></p> <p>Allow the water to run for 3 to 5 minutes prior to taking the sample to flush the water lines. Decrease the water flow to the diameter of a pencil to reduce splashing.</p> <p>Do not rinse the sample bottle prior to sampling; it contains a chemical preservative. Fill the sample 1/3 full. <b>Caution: Hydrochloric Acid (HCl) is corrosive</b> Carefully add three drops of the 1:1 Hydrochloric acid. Fill the sample vial completely full (zero headspace) without overflowing and flushing out the preservatives. Be sure the cap is tightened to prevent leakage during shipment to the laboratory.</p> <p>Sample must be collected in duplicate.</p> <p><b>RECORD BOTTLE NUMBER ON FORM</b></p>
(TTHM) Total Trihalomethanes only	Cool 6°C, (0.003 g) Sodium Thiosulfate	14 days	40 ml	Glass with Teflon Septa 	<p><b>TTHM only Sampling Instructions</b></p> <p>14 Day Holding Time Overnight Delivery Shipped Cold &lt;6°C</p> <p>Use the cold-water faucet. <b>Remove any faucet attachments and the aerator prior to sampling. Also, avoid rubber hose attachments.</b></p> <p>Allow the water to run for 3 to 5 minutes prior to taking the sample to flush the water lines. Decrease the water flow to the diameter of a pencil to reduce splashing.</p> <p>Do not rinse the sample bottle prior to sampling; it contains a chemical preservative. Fill the sample vial completely full (zero headspace) without overflowing and flushing out the preservatives. Be sure the cap is tightened to prevent leakage during shipment to the laboratory.</p> <p>Sample must be collected in duplicate.</p> <p><b>RECORD BOTTLE NUMBER ON FORM</b></p>