

Office of Laboratory Services 167 11<sup>th</sup> Avenue South Charleston, WV 25303 Phone: (304) 558-3530 | FAX: (304) 558-2006

# Instructions for the Supply Order Form for Diagnostic Immunology Test Kits

## NOTE: This is NOT a supply order form

Supplies are provided at no charge to facilities that submit the specimens to the Office of Laboratory Services.

#### **FORM INSTRUCTIONS:**

FROM:							
MAILING ADDRESS:							
CITY:			STATE:	ZIP:			
PHONE: EN		EMAIL:					
NAME:	TITLE:			DATE:			

#### Please print legibly.

- From: The complete name of your facility. Do not use the initials of your site.
- Mailing Address: Include the full address of your facility.
- Supplies cannot be shipped to a Post Office Box.
- Phone: Include area code.
- **Name and Title**: Include the name of the person responsible for ordering supplies from your facility. This will allow the lab to contact them directly if there is a problem with the order.
- Date: Include the date that the order was placed.

#### SYPHILIS, RUBELLA, HIV OR HEPATITIS TESTING SUPPLIES:

Syphilis, Rubella, HIV or Hepatitis Testing Supplies	Number Ordered	Number Sent
Blood Tubes [8.5ml BD Vacutainer Tubes]		
Needle Hubs [Eclipse Needle 21G w/ holder]		

#### • Blood Tubes:

• Write the requested number of blood tubes in the "Number Ordered" portion of the form.

Instructions for the Supply Order Form for Diagnostic Immunology Test Kits - 1

- OLS provides 8.5mL BD Vacutainer Tubes also called SST tubes or tiger type tubes.
- NOTE: OLS requires one (1) blood tube per test when requesting Syphilis, Rubella or HIV test and two
  (2) blood tubes when requesting Hepatitis B and C.

## • Needle Hubs:

- Write the requested number of needle hubs in the "Number Ordered" portion of the form.
- OLS provides the Eclipse needle 21gage with holder. This allows clinicians to draw multiple blood tubes but only stick the patient once.
- Facilities should order one (1) hub for every two (2) blood tubes ordered.

## CT/NG (URINE) SUPPLIES:

- Urine Collection Device:
  - Write the number of urine collection devices requested in the "Number Ordered" portion of the form.
  - OLS provides the Hologic Aptima<sup>®</sup> collection device. The collection device is pre-packaged by Hologic in an easy open packet containing a collection device and plastic transfer pipette.

## • Urine Cups:

• OLS provides 4.5-ounce specimen cups with screw top lids; one for each collection device ordered.

### **ADDITIONAL SUPPLIES:**

Additional Supplies	Number Ordered	Number Sent
Specimens Submission Forms		1 Hepatitis Form
Newest form available online: www.dhhr.wv.org/ols		1 Blood/Urine Form
Shipping canisters with 2 absorbent pads (holds 8 tubes)		

#### • Specimen Submission Forms:

- Hepatitis Form:
  - The Hepatitis Form is specifically designed for ordering Hepatitis A, B, or C.
  - One (1) copy of the form and a set of instructions with more information will be included with each order.
- o Blood/Urine Form:
  - The Blood/Urine form is designed for ordering CT/NG, HIV, Rubella, Syphilis, or HIV Confirmation for those sites enrolled in the HIV Rapid Program.

- One (1) copy of the form and a set of the instructions with more information, will be included with each order.

## • Shipping Canisters:

- Write the number of shipping canisters requested in the "Number Ordered" portion of the form.
- Packaging provided by OLS meets all current Postal and Department of Transportation (DOT) shipping regulations.
- Each canister assembly comes with two (2) absorbent pads and holds up to eight (8) blood or urine tubes.
- At this time, the facility is responsible for cost of shipping samples to the laboratory.